


Transmitters for R-SW Handpiece »FALCON«



Radial Transmitters – Enthesopathies/Trigger Points

A6
»AkuSt«



- Ø 6 mm
- Shock wave acupuncture
- Penetration depth: 0 – 20 mm
- Intensity:

F15
»White Soft«



- Ø 15 mm
- Superficial pain regions, muscles of mastication, cervical spine
- Penetration depth: 0 – 30 mm
- Intensity:

T10
»Finger«



- Ø 10 mm
- Trigger points, deep-lying enthesopathy
- Penetration depth: 0 – 40 mm
- Intensity:

R15
»Basic Energy«



- Ø 15 mm
- Any type of enthesopathy
- Penetration depth: 0 – 35 mm
- Intensity:

D115
»Deep Impact®«



- Ø 15 mm
- Deep-lying pathologies, chronic conditions, high-energy application
- Penetration depth: 0 – 50 mm
- Intensity:

C15
»CERAmA-x®«



- Ø 15 mm
- Any type of enthesopathy
- Penetration depth: 0 – 30 mm
- Intensity:

D-ACTOR® Transmitters – Myofascial Treatment

D20-S
»Basic Oscillator«



- Ø 20 mm
- Muscle tissue and connective tissue
- Penetration depth: 0 – 50 mm
- Intensity:

D20-T
»Golden Oscillator«



- Ø 20 mm
- For professional athletes, high-energy application
- Penetration depth: 0 – 50 mm
- Intensity:

D35-S
»Large Oscillator«



- Ø 35 mm
- Muscle tissue and connective tissue, large muscles
- Penetration depth: 0 – 60 mm
- Intensity:

Special Transmitters – Enthesopathies/Trigger Points

Ro40
»Beam«




- Ø 15 mm
- Any type of enthesopathy
- Penetration depth: 0 – 35 mm
- Intensity:

B15
»White Tip«



- Ø 15 mm
- Any type of enthesopathy
- Penetration depth: 0 – 30 mm
- Intensity:

ATLAS
»Soft Tip«
(after Di Maio/Swart)



- Ø 15 mm
- Any type of enthesopathy
- Penetration depth: 0 – 10 mm
- Intensity:

Transmitters for R-SW Handpiece »FALCON«



PERI-ACTOR® Transmitters (according to Swart/Di Maio) – Fascial Treatment

PERI-ACTOR®
»scraper«



■ Treatment of all fascia of the musculoskeletal system

PERI-ACTOR®
»knuckle«



■ Treatment of all medium-sized and large fascia

PERI-ACTOR®
»sphere«



■ Spot treatment of fascia

PERI-ACTOR®
»scoop«



■ Treatment at or under bone edges



SPINE-ACTOR® Transmitters – Treatment of the Paravertebral Muscles

SPINE-ACTOR® 1



■ Span: 20 mm

SPINE-ACTOR® 2



■ Span: 27.5 mm

SPINE-ACTOR® 3



■ Span: 35 mm

The SPINE-ACTOR® transmitters have been developed to treat the paravertebral muscles. They come in three sizes allowing treatment in all sections of the spinal column muscles (lumbar, thoracic and cervical).

